



Developed in the UK  
by Growth Technology

**Growth  
Technology**

**0.00059 - 0.00004 - 0.00005**

**GUARANTEED ANALYSIS:**

Total Nitrogen (N) ..... 0.00059%  
0.00059% Nitrate Nitrogen

Available Phosphate ( $P_2O_5$ ) . 0.00004%

Soluble Potash ( $K_2O$ ) ..... 0.00005%

*Derived from: calcium nitrate  
and dipotassium phosphate*

*Information regarding the  
contents and levels of metals  
in this product is available  
on the internet at: [http://  
www.regulatory-info-hd.com](http://www.regulatory-info-hd.com)*

04US

F-1507

**Net Weight**  
**100 g / 3.53 oz.**



5711 Enterprise Drive  
Lansing, MI 48911

**Growth Technology**

# CLONEX<sup>®</sup> MIST

**Root  
Promoter™**

**100 ml / 3.38 fl oz.**



## **Instructions for Use:**

If possible spray the tips of the mother plant with **Clonex Mist** once or twice in the two weeks before taking cuttings.

We recommend **Root Riot** cubes but **Rockwool** cubes are also useful.

**1/** Take cuttings by making a clean diagonal cut just below a leaf node.

**2/** Dip into **Clonex Rooting Compound** immediately and insert cutting into a rooting cube.

**3/** Spray new cuttings with **Clonex Mist** and place them into a propagator.

**4/** After 2–3 days mist the cuttings again.

**5/** Repeat misting every two days, or whenever cuttings look dry.

**[www.hydrodynamicsintl.com](http://www.hydrodynamicsintl.com)**